



All ceramic crown

We are Chinese dental lab. We provide custom-made dental service for dental institutions at home and abroad, and provide oem service for other dental lab. All ceramic crown is one of our main products. Made from materials certified by FDA and the European Union. We supply full contour zirconia, layered zirconia porcelain, multilevel full contour zirconia and emax. Crowns, bridge, veneer, inlay and onlay will be the best choice for you.

All Ceramic Crown Outsourcing Dental Lab in China provides lab to lab service



Shenzhen Wanmei Dental Laboratory Co.,LTD.

Tel: +86-755-29760054

E-mail: info@wmdentallab.com



As for aggressive selling prices, we believe that you will be searching far and wide for anything that can beat us. We are able to state with absolute certainty that for such good quality at such costs we've been the lowest around for Low price for China All ceramic crown Dental Lab, "Quality", "honesty" and "service" is our principle. Our loyalty and commitments remain respectfully at your provider. Make contact with Us Today For more data, make contact with us now.

Low price for China Orthodontic Dental Laboratory, Orthodontic Dental Labs, We integrate design, manufacture and export together with more than 100 skillful workers, strict quality controlling system and experienced technology. We keep long term business relationships with wholesaler and distributors form more than 50 countries, such as USA, UK, Canada, Europe and Africa etc.

Just as their name implies, all-ceramic dental crowns are restorations whose full thickness is made entirely of a glass-like substance, such as porcelain.

BENEFITS:

- I Metal Free and Durability
- I Virtually Chip-proof (up to 1200 MPa)
- I Milling by 5 axis machine.
- I Five different zirconium restoration types we provide.
- I Designed and milled using CAD / CAM technology.
- I Life-time warranty for full-contour zirconia crown.
- I Perfect biocompatibility, No stimulation to gum tissue.

Shenzhen Wanmei Dental Laboratory Co.,LTD.

Tel: +86-755-29760054

E-mail: info@wmdentallab.com



Turnaround/Shipping

All cases are 4days in the lab only! Send back by DHL、UPS or FEDEX!



In voice

The invoice will send to you by email after cases send back



Statement

Monthly statements will be sent out on the last day of each month



Payment Terms

Next 15 days from statement date, or no more than \$5000



Warranty

works are redo/repair within 2 years due to defective material or workmanship.
*does not cover breakage due to negligence



Methods of Payment

Bank transfer/Paypal



Material List

Type	Name of material	CE No.	Lot No.	Manufacturer	Composition
Non-precious Alloy	Argeloy N.P.Supreme	CE0086	14034132	Argen(USA)	Co60.%,Cr27%,W5%,Mo6%
Precious Alloy	Argedent 74	CE0086	718846622	Argen(USA)	Au74%,In2.5%,Pt4.0%,Pd7.6%
	Argedent 52	CE0086	71977557	Argen(USA)	Au52.5%,Pd26.9%,Ag16%,In2.5%
Semi-precious	Argistar 43+	CE0086	7545643	Argen(USA)	Au2%,Pd42.9%,Ag46%,In8%
Metal Framework	Argeloy N.P.Partial	CE0086	13869812	Argen(USA)	Co64%,Cr28%,Mo6%
	Vitalium 2000	CE01663	N001830	Dentsply(USA)	Co63.8%,Cr28%,Mo5.9%
Acrylic Resin	TCS Unbreakable	CE0470	73021314	TCS(USA)	Acrylic Resin
	Vertex	CE120	XT016907	Castaperess	Acrylic Resin
	Lucitone 199	CE PR017600	131114	Dentsply(USA)	Acrylic Resin
All Ceramic	E max Press 3	CE0123	P19900	Ivoclar(Liechtenstein)	SiO279%,LiO15%,K2O5%,P2O5%,ZrO2%,ZrO3%
	Zirconia	HD601143840001	W1802260NG-1	Aidite(China)	ZrO2,Y2O3,AL2O3
	Full Contour Zirconia	HD601143840001	W171023QGC-1	Aidite(China)	ZrO2,Y2O3,AL2O3
Ceramic Powder	E max Ceram 3	CE0123	P62370	Ivoclar(Liechtenstein)	SiO265%,Al2O38%,Na2O6%,ZrO3%
	Vita Ceram	CE0124	LOT45190	Vita(Germany)	SiO265%,Al2O315%