



### **Bionator Appliance**

The following is about Bionator Appliance related, I hope to help you better understand Bionator Appliance.

**Bionator Appliance outsourcing dental lab in China with lab to lab service**





Having a positive and progressive attitude to customer's fascination, our organization constantly improves our solution high-quality to fulfill the requirements of shoppers and further focuses on safety, reliability, environmental prerequisites, and innovation of Chinese Professional Bionator Appliance From China Dental Laboratory, We have been also the appointed OEM factory for several worlds' famous merchandise brands. Welcome to speak to us for more negotiation and cooperation.

Chinese Professional China Bionator Appliance, Pendulum Orthodontic Appliance, So We also continuously function. we, focus on high quality, and are conscious of the importance of environmental protection, most of the merchandise are pollution-free, environmentally friendly solutions, reuse on the solution. We've Updated our catalog, which introduces our organization. n detail and covers the primary products we provide at present, You may also visit our web-site, which involves our most recent product line. We look forward to reactivating our company connection.

The Bionator I (to open bite) is designed to aid in correction of Class II malocclusions by maintaining the mandible in an advanced position and guiding eruption of the posterior teeth. In order to prevent the super eruption of the anterior teeth, the mandibular anteriors are covered with an acrylic cap which comes in contact with the maxillary anteriors. An expansion screw can be added to maintain a tight fit of the appliance.

Shenzhen Wanmei Dental Laboratory Co.,LTD.

Tel: +86-755-29760054

E-mail: [info@wmdentallab.com](mailto:info@wmdentallab.com)

The Bionator II (to close bite) is designed to aid in correction of Class II malocclusions and the reduction of an anterior open bite. The latter is achieved by keeping the acrylic coverage over the occlusal surfaces of the posterior teeth. An expansion screw can be added to maintain a tight fit of the appliance.

Hot Tags: Bionator Appliance outsourcing dental lab in China, Bionator Appliance manufacturer, Bionator Appliance lab to lab service, Dental crowns, Zirconia crown, PFM crown, Dental implants dental center